


Prosigna®

Breast Cancer Assay

2nd generation test
that can help physicians more comprehensively
inform treatment decisions.



**More
Accurate
Prognostic
Score**

**Intrinsic
Molecular
Subtype
Identified**

**Clinicopathological
Factors
Integrated**



Prosigna®

Breast Cancer Assay

Prosigna® helps to accurately distinguish high-, intermediate-, and low-risk groups¹ when compared to the 1st generation assays and can better inform treatment decisions.¹

2nd GENERATION ASSAY Prosigna®

- 50 genes
- Gene expression + clinical factors (age, node status, tumor size, histologic grade)
- Results can be available in as little as 24 hours


10-YEAR PROGNOSIS

1st GENERATION ASSAY

- 21 genes
- Clinical and pathological factors not included in risk calculation
- Samples need to be sent out

5 to 9-YEAR PROGNOSIS


	Number of Women	Number of DR	10-year DR Risk (%)
Prosigna ROR low and Oncotype Dx RS low	104	2	2.1%
Prosigna ROR low, but Oncotype Dx RS intermediate or high	261	7	2.9%
Prosigna ROR intermediate or high, but Oncotype Dx RS low	62	10	18.2%
Prosigna ROR high and Oncotype Dx RS high	55	14	19.5%



With 1st generation assays, nearly 20% of patients can be wrongly assigned to a low-risk group and not receive life-saving chemotherapy.¹

Clinically Validated Algorithm Generates a Prosigna® Score Unique for Each Patient.

Prosigna® measures the expression of 50 different genes to identify subtype and report a Risk of Recurrence Score (ROR), which is used to assign the patient to a predefined risk group. These results presented in the Patient Report are derived from a proprietary algorithm based on the PAM50 gene signature, intrinsic subtype, and clinical variables including tumor size and nodal status.




How does the Prosigna® Assay work?

1

Extract RNA from FFPE* tumour sample

with the Veracyte FFPE RNA Extraction Kit



2

Test RNA using Prosigna® Assay


on the nCounter Analysis System in your local lab in as little as 24 hours



3

Capture expression profile

of patient's tumour for tailored, more accurate and better informed treatment decisions




4

Patient Report

to enable fast and precise treatment decisions:

- **Intrinsic Molecular Subtypes identified**
- **Risk of Recurrence (ROR) within 10 years**
- **Clinicopathological factors** like tumour size and nodal status incorporated in the score



*FFPE: Formalin Fixed Paraffin Embedded Tissue.


Prosigna®

Breast Cancer Assay

More **comprehensive information about your patient's tumour** helps provide fast, precise and tailored **treatment decisions**.




Provides more accurate prognosis to help improve breast cancer treatment and outcomes¹.



Combines tumour gene expression and clinicopathological factors in a single 10-year Risk of distant Recurrence (ROR) score^{2,3} for informed treatment decisions.



The only breast cancer prognostic test identifying the four PAM50 intrinsic molecular subtypes for tailored treatments.



Easy access to local testing with faster turnaround times. Results can be available in as little as 24 hours.

References:

1. Sestak I, Buus R, Cuzick J, et al. JAMA Oncol. 2018; 4(4):545-553. 2. Kos et al. Breast Cancer Res, 2014, 16:103. 3. Alexandre et al. Cancer Manag Res, 2019; 11: 10353-10373.